

Coagulation Testing Update

Effective January 3rd, 2023

Effective January 3, 2023, Wisconsin Diagnostic Laboratories (WDL) will no longer be performing PF4 IgG for heparin antibodies in-house. PF4 IgG testing is currently being performed at Versiti. As of January 3, 2023 the in-house PF4 IgG assay will be no longer be orderable in Epic and will be replaced by the Versiti PF4 IgG assay.

Discontinued Tests:

Test	Mnemonic	Orderable	CPT code
Platelet Factor 4 for Heparin Antibodies IGG	PF4	LAB10502	86022

Versiti send out test replacing the discontinued orderable:

Test	Mnemonic	Orderable	CPT code
HIT Evaluation (PF4 IgG With Reflex to SRA)	HIT / PF4	LAB11149	86022

Key points for the Versiti test:

- Turnaround time is 1-2 days but testing is not performed on weekends or holidays
- If the preliminary result is positive, the confirmatory SRA test will be done
- Serum is the required specimen type – 3 red top tubes should be collected and 3 frozen aliquots should be sent

For questions, please contact Wisconsin Diagnostic Laboratories Coagulation department:

- Alexandra Harrington, MD, MLS(ASCP)^{CM}, Medical Director, Hematology and Coagulation, 414.805.6703, aharrington@mcw.edu
- Ashley Cunningham, MD, Associate Medical Director, Coagulation, 414.805.6943, acunningham@mcw.edu
- Chuck Schur, Manager - Hematology, Coagulation and Flow Cytometry, 414.805.5159, CSchur@wisconsindiagnostic.com
- Kenneth Kohlhapp, Technical Specialist – Coagulation, 414.805.7649, kkohlhapp@wisconsindiagnostic.com

